

EDISON 사이언스 앱을 사용한 비구획분석과 생물학적동등성 분석의 통합

서울아산병원
임상약리학과
전공의 한성필



Introduction: Bioequivalence


- Definition
 - When administered at the **same molar dose** under similar conditions in an appropriately designed study
 - **Absence of a significant difference**
 - At the site of drug action
 - The **rate and extent** to which the **active ingredient or active moiety** in pharmaceutical equivalents or pharmaceutical alternatives

Introduction: Generic Drug

- Drug product which is identical or bioequivalent to Reference drug in:
 - Active ingredient (s)
 - Route of administration
 - Dosage form
 - Strength
 - Indications
 - Safety
- May have different:
 - Inactive ingredients
 - Colour
 - Shape
- Almost half of drugs in market have Generics

Introduction: Bioequivalence test

- 생물학적동등성시험(bioequivalence test)는 기존 의약품의 특허가 만료된 후, 해당 의약품을 동일하게 개발하여 판매하고자 할 때 수행하는 임상시험.



(Source: 식약처 가이드라인)

Introduction: Bioequivalence test

- Study Design
 - Good experimental design, enhances the power of the study
 - Depends on: question to be answered, nature of reference drug/ dosage form, benefit-risk ratio
- As far as possible, the study should be of crossover design & suitably randomized
- Ideal design:
 - Randomized two-period, two-sequence, crossover design with adequate washout period


Introduction: Parameters to be measured

- Pharmacokinetic Parameters measured are:

- C_{\max}
- AUC_{0-t}
- $AUC_{0-\infty}$
- (T_{\max})

- For steady state studies:

- AUC_{0-t}
- C_{\max}
- C_{\min}
- Degree of fluctuation





Introduction: Assessment

- PK Parameter
 - log-normal distribution
- Assessment
 - Log-transformed C_{\max}
 - Log-transformed AUC_{0-t}
 - geometric mean ratio – 90% confidence Interval
 - $0.80 \sim 1.25$


서론

- 각각의 소프트웨어가 사용하기 어렵고 과정이 복잡하며 큰 비용이 소모. 특히 약동학 분석 초보자나 통계학 비전공자가 사용하기 어려움. 본 저자들은 EDISON 사이언스 앱을 사용하여 이 과정을 연속성을 지닌 한 과정으로 통합.
- 본 연구가 제시하는 방법을 통해 쉽고 정확하고 비용이 들지 않는 빠른 비구획분석과 생물학적 동등성 분석이 가능하다는 것을 보이고자 함.

Dataset: 2x2 BE test




EDISON Science Apps




PK 1 Comp Oral

My EDISON > pk1coma



PK 1 Comp IV

My EDISON > pk1cima




NonCompartmentEdison: Noncompartmental Analysis for Pharmacokinetic Parameters using...



생물학적동등성 시험 자료 분석 및 시뮬레이션

My EDISON > edisonBE

EDISON Science Apps




PK 1 Co

My EDISON > pk1

PK 1 Comp IV

My EDISON > pk1cima



NonCompartmentEdison
using...



생물학적동등성 시험

My EDISON > edisonBE

nSubj : 10

CL : 30

V : 100

Ka : 2

DH1 : 0, 100000

Time : 0, 0.25, 0.5, 1, 2, 4, 5, 7, 9, 12, 24

PropE : 0.1

AddE : 0.1

LLoQ : 0

Jitter : 1


FullCov : 0.04, 0.03, 0.03, 0.04

concLog : Linear ▼


Pharmacokinetic Parameters

Simulated Data

Treatment — R — T



EDISON Science Apps




PK 1 Comp Oral

My EDISON > pk1coma



PK 1 Comp IV

My EDISON > pk1cima



NonCompartmentEdison: Noncompartmental Analysis for Pharmacokinetic Parameters using...



생물학적동등성 시험 자료 분석 및 시뮬레이션

My EDISON > edisonBE


입력 자료의 형태

			SUBJ	GRP	PRD	TRT	nTIME	TIME	CONC
2x2 BE trial	Period 1	Period 2	1	RT	1	R	0	0	0
N=12			1	RT	1	R	0.25	0.26	511.3
Sequence 1 (BA)	Y_{11}	Y_{12}	1	RT	1	R	0.5	0.46	678.79
BA is RT			1	RT	1	R
Sequence 2 (AB)	Y_{21}	Y_{22}	1	RT	2	T	0	0	0
AB is TR			1	RT	2	T	0.25	0.25	487.62
			1	RT	2	T	0.5	0.48	769.6
		
			5	TR	1	T	0	0	0
			5	TR	1	T	0.25	0.23	382.79
			5	TR	1	T	0.5	0.45	477.03
			5	TR	1	T
			5	TR	2	R	0	0	0
			5	TR	2	R	0.25	0.28	596.98
			5	TR	2	R	0.5	0.47	832.76
			5	TR	2	R
		

결과 자료의 형태


SUBJ	GRP	PRD	TRT	AUClast	Cmax	Tmax
1	RT	1	R	5018.927	1043.13	1.04
1	RT	2	T	6737.507	894.21	1.03
2	TR	1	T	4373.97	447.26	1.01
2	TR	2	R	6164.276	783.92	1.98
4	TR	1	T	5592.993	824.42	1.97
4	TR	2	R	5958.16	646.31	0.97
5	TR	1	T	3902.59	803.7	0.8
5	TR	2	R	4620.156	955.3	0.74

EDISON Science Apps




PK 1 Comp Oral

My EDISON > pk1coma



PK 1 Comp IV

My EDISON > pk1cima



NonCompartmentEdison: Noncompartmental Analysis for Pharmacokinetic Parameters using...



생물학적동등성 시험 자료 분석 및 시뮬레이션

My EDISON > edisonBE

- Comparison of 90% confidence interval for the ratio of the geometric means of (A) AUC_{0-t} and (B) C_{max}

(A)

Analysis	Lower Limit	Point Estimate	Upper Limit
EDISON Science App	0.88944	0.95408	1.02341
SAS: PROC GLM	0.88944	0.95408	1.02341
SAS: PROC MIXED	0.88944	0.95408	1.02341

(B)

Analysis	Lower Limit	Point Estimate	Upper Limit
EDISON Science App	0.90136	0.97984	1.06515
SAS: PROC GLM	0.90136	0.97984	1.06515
SAS: PROC MIXED	0.90136	0.97984	1.06515

결과

1. ANOVA Table

```
> print(test2x2(NCAREsult4BE, "AUClast"), na.print="")
$`Analysis of Variance (log scale)`
      SS DF      MS      F      p
SUBJECT 2.875497e+00 32 8.985928e-02 3.183942248 0.0008742828
GROUP   1.024607e-01  1 1.024607e-01 1.145416548 0.2927731856
SUBJECT(GROUP) 2.773036e+00 31 8.945279e-02 3.169539016 0.0009544080
PERIOD   3.027399e-05  1 3.027399e-05 0.001072684 0.9740824428
DRUG     3.643467e-02  1 3.643467e-02 1.290972690 0.2645764201
ERROR    8.749021e-01 31 2.822265e-02
TOTAL    3.786834e+00 65
```

2. Variability

```
$`Between and Within Subject Variability`
      Between Subject Within Subject
Variance Estimate      0.03061507      0.02822265
Coefficient of Variation, CV(%) 17.63193968 16.91883011
```

3. LSM

```
$`Least Square Means (geometric mean)`
      Reference Drug Test Drug
Geometric Means      5092.098 4858.245
```

4. 90% CI of GMR

```
$`90% Confidence Interval of Geometric Mean Ratio (T/R)`
      Lower Limit Point Estimate Upper Limit
90% CI for Ratio 0.889436      0.9540753 1.023412
```

5. Sample size

```
$`Sample Size`
      True Ratio=1 True Ratio=Point Estimate
80% Power Sample Size      6      7
```

Conclusion

- 본 연구는 EDISON 사이언스 앱을 사용해 쉽고 정확하고 비용이 들지 않는 빠른 비구획분석과 생물학적동등성 분석법을 제시하였다.
- 현재는 이러한 분석을 위해서 여러 상용 소프트웨어를 필요로 하는 복잡한 단계를 거쳐야 한다.
 - 따라서 분석 시간이 오래 걸리고 많은 비용이 소모되었다.
 - 본 저자들은 EDISON 사이언스 앱을 사용하여 이 두 과정을 통합하였고, 이로서 농도-시간 자료로부터 비구획분석과 생물학적동등성 통계 분석까지 연속적으로 가능케 하였다.
- 본 분석에서 학계와 산업계에서 가장 정확하게 생물학적 동등성을 평가한다고 말하는 SAS통계 패키지와의 비교를 통해 EDISON 사이언스 앱을 사용한 생물학적 동등성 분석이 정확한 값을 제시함을 검증할 수 있었다.